










Antikoagülan Tüpler (EDTA'lı, Heparinli, Sitratlı), Plastik PP



Vakumsuz

Hacim	3 mL	5 mL	5 mL	5 mL	5 mL	5 mL	5 mL	5 mL	8 mL	
Çap x Boy (mm)	12 x 51	12 x 75	16 x 60	12 x 92	12 x 75	16 x 60	12 x 92	16 x 75	16 x 100	
	KF + Na ₂ EDTA		L115317 2.5 mL kan için							
	K2EDTA	L113337 2.5 mL kan için	L110437 2.5 mL kan için	L114337 2.5 mL kan için	L115337 2.5 mL kan için	L116337 5 mL kan için	L117337 5 mL kan için	L118337 5 mL kan için		
	K3EDTA	L113397 2.5 mL kan için	L110487 2.5 mL kan için	L114397 2.5 mL kan için	L115397 2.5 mL kan için	L111497 5 mL kan için	L116397 5 mL kan için	L117397 5 mL kan için	L117497 5 mL kan için	L118397 5 mL kan için
	Lityum Heparin	L113357 2.5 mL kan için	L110457 2.5 mL kan için	L114357 2.5 mL kan için	L115357 2.5 mL kan için	L116357 5 mL kan için	L117357 5 mL kan için	L118357 10 mL kan için		
	Sodyum Sitrata	L114867 0.4 mL		L115867 0.4 mL	L114367 0.5 mL	L115367 0.5 mL	L118367 0.5 mL			
	Sodyum Sitrata	L113577 0.25 mL	L114577 0.25 mL	L115577 0.25 mL						
	Sodyum Sitrata	L113877 0.4 mL	L114877 0.4 mL	L115877 0.4 mL						
	Sodyum Sitrata	L113477 0.5 mL	L114477 0.5 mL	L115477 0.5 mL	L114377 1 mL	L115377 1 mL	L118377 1 mL			