

## Cam Pamuğu

| Katalog No | Ambalaj |
|------------|---------|
| C490001    |         |



## Tıpa, Şilif için

### Cam

Borosilikat 3.3 camdan, tutma yeri hegzagonal.

| Katalog No | Şilif (NS) | Ambalaj / Adet |
|------------|------------|----------------|
| C482000    | 7 / 16     |                |
| C482001    | 10 / 19    |                |
| C482002    | 12 / 21    |                |
| C482003    | 14 / 23    | 10             |
| C482004    | 19 / 26    |                |
| C482005    | 24 / 29    |                |
| C482006    | 29 / 32    |                |



## Plastik PE, Konik

| Katalog No | Şilif (NS) | Ambalaj / Adet |
|------------|------------|----------------|
| L291080    | 7 / 16     |                |
| L291081    | 10 / 19    |                |
| L291082    | 12 / 21    | 50             |
| L291083    | 14 / 23    |                |
| L291084    | 19 / 26    |                |
| L291085    | 24 / 29    |                |
| L291086    | 29 / 32    |                |
| L291087    | 34 / 35    | 10             |
| L291088    | 45 / 40    |                |
| L291089    | 60 / 50    | 1              |



## Teflon PTFE, Konik

| Katalog No | Şilif (NS) | Uzunluk (mm) |
|------------|------------|--------------|
| L295640    | 10 / 19    | 31           |
| L295641    | 14 / 23    | 34           |
| L295642    | 19 / 26    | 38           |
| L295643    | 24 / 29    | 44           |
| L295644    | 29 / 32    | 47           |
| L295645    | 34 / 35    | 51           |

